




Innovate Dentistry with Eighteeth

www.eighteeth.com

 Changzhou Sifary Medical Technology Co., Ltd
Add: NO.99, Qingyang Road, Xuejia County, Xinbei
District, Changzhou City, 213000 Jiangsu, China
Email: ivy@sifary.com
Web: www.eighteeth.com

DISTRIBUTED BY



CATALOG

1

Endodontics & Restorative

FindPex Wireless Apex Locator	01-02
E-Pex Apex Locator	03-04
AirPex Wireless Apex Locator	05-06
E-STYLE & E-FLEX Series NiTi Endo File	07-08
E-xtreme Cordless Endo Motor	09-10
E-CONNECT Cordless Endo Motor	11-12
E-connect S Endomotor With Built-in Apex Locator	13-14
E-connect S+ Endomotor with Built-in Apex Locator	15-16
E-value Cordless Endo Motor	17-18
Ultra X Cordless Ultrasonic Activator	19-20
Ultra Mint & UltraMint Pro All-round Ultrasonic Scaler	21-26
Fast-Pack pro & Fast-Fill Cordless Obturation System	27-34
Gutta Percha Points	35
Absorbent Paper Points	36
CuringPen 4 LEDs Curing Light	37-38
CuringPen-E 3 LEDs Curing Light	39-40

2

X-Ray Imaging Systems

HyperLight Portable X-ray Unit	41-44
HyperLight-G Portable X-ray Unit	45-46
NanoPix Intraoral X-Ray Sensor	47-50
NanoPix-E 1.3 Intraoral X-Ray Sensor	51-52
FinScan F350 Oral and Maxillofacial CBCT Equipment	53-56

3

Digital Impressions

Helios 500 Intraoral Scanner	57-62
Helios 600 Intraoral Scanner	63-68

4

Optics

AcuVision X Surgical Microscope	69-76
Brilliance Dental Loupes	77-78
Brilliance Business Dental Loupes	79-80

CATALOG

Brilliance 48° 81-84

Dental Loupes

Softouch 85-86

Surgical Headlight

Wireless Z 87-88

Wireless Headlight

5 Implant & Surgery Systems

MotorSurg 89-92

Dental Implant Unit

E-FLOW 93-94

Dental Anesthesia Device


MotorTurbo & E-Asp 1 95-98


Electric Motor and 1:5 Handpiece

6 Large Equipment

E-Sanit 99-100

Steam Sterilizer





FindPex

Apex Locator Using ASW Technology



Main Technical Specifications:

Dimensions	156mm × 106mm × 114mm (package)
Weight	247g
Screen	3.5" color TFT
Battery	Lithium-ion battery, 3.7V/800mAh
Adapter	Input: AC 100-240V, 50Hz/60Hz, 0.2A Output: DC 5V/1A
Frequency	50Hz/60Hz, ±10%
Power rating	0.2A MAX
Degree of protection	IPX 0
Electrical safety class	Class II




More Accurate & Anti-Interference

- Exclusive **Adaptive Sine Wave (ASW)** technology & next-gen multi-frequency technology enhance the measurement accuracy.
- High-speed AI chip allows real-time response
- Adaptive to various root canal conditions


3 Adjustable Measurement Sensitivity Levels

Tailored to meet the demands of different users



Personalized Customization


- Personalized background & Sound recognition
- Advanced settings include time for auto-shut down & auto-return to the home page, language and volume
- Adjustable screen brightness fits different visual needs & saves energy



Light and Dark Theme Interface

Ultimate Experience of Visual & Tactile Sensations


- Screen visible from any angle: 3.5-inch TFT color screen and IPS technology



- Upgraded touch buttons

Intuitive Operation


- Synchronous endo file position demonstration
- Distinguishes Oral & Non-Oral Environments: 'Connected' indicator for non-oral settings, numerical display for oral environments



Long-lasting Durability

- Durable Accessories: 316L stainless steel lip hook & biomimetic-silica-wrapped file clip





E-PEX

Smart Apex Locator



Integrate Endomotor with Apex Locator (Connect Function)

- Auto Start & Stop
- Apical Reverse
- Apical Slow Down




Smart Accurate and Reliable Apex Locator

- Advanced Multi - frequency Technology
- High Precision in Wet & Dry Canals
- Large 3.5" LCD Display
- Automatic Calibration
- Compact Size & Stable Design

Screen Display & Its Meaning

	Interface	Interface	Interface	Interface
Display				
Value	1 2	0 9	0 0	- 3
Meaning	The corresponding scale value is 2, indicating the distance from the apex is about 2mm	It indicate that the distance from the apex is about 1.5mm	It indicates that the endo file reaches the apex	It indicates that the root canal file has exceeded the apex
Alarm	Didi alarm sound for a long time interval	Didi alarm sound interval becomes shorter	Uninterrupted long beep	Very rapid Didi alarm sound




AirPex

Wireless Apex Locator

International Patent

No.:PCT/CN2019/110429




Feel the Charm of Wireless Charging

Lay it down , Charge it up
Charging is as simple as placing
the AirPex on the charging base



True Wireless Apex Locator

Get rid of the wire to the fullest
Use without magnet wire
Accuracy through wirelessness





E-STYLE & E-FLEX Series

NiTi Endo File



- View Details -



E-FLEX RT



E-FLEX S



E-STYLE BLUE

E-FLEX REC



E-FLEX ONE



E-STYLE GOLD

E-FLEX GOLD



E-FLEX BLUE



E-FLEX MINI



E-FLEX EDGE





E-xtreme

Cordless Endo Motor

Smaller | Lighter | Quieter than ever



360° contra-angle rotation



Slim and compact design for maximum comfort and mobility



Miniature contra-angle for better visibility & easier posterior access

– Miniature contra-angle Super Mini head & thin neck



– Normal head & neck




Eighteenth innovative motor, cutting down the noise from inside, allowing you to experience the tranquility.



Less than 100 grams, reducing hand fatigue during treatment, providing the ultimate E-xtreme experience.



Reciprocation function



Built-in file system


Specifications

Model	: E-xtreme	Power supply	: 100 V-240 V
Weight(grams)	: 99 g	Frequency	: 50 Hz / 60 Hz
Torque range	: 0.5 N.cm-4.0 N.cm	Output	: DC 5 V/ 1 A



E-CONNECT

Cordless Endo Motor



- Torque Calibration
- Powerful 1500 mAh Lithium Battery
- Digital OLED Display
- Integration Option
- Low Noise Motor
- Screen Cast / Mirroring
- Ambidextrous * Display
- *User can rotate the Endomotor Screen towards left or right
- *In connect function Endomotor mirrors Apex Locator's screen

Integrate Endomotor with Apex Locator (Connect Function)

- Auto Start & Stop
- Apical Reverse
- Apical Slow Down




Five Reciprocation Modes (M0)

- Reverse Reciprocation-One Mode**
Counter clockwise & clockwise movement
Fwd 30° - Rev 150°
For reciprocating files
- Forward Reciprocation-Four Modes**
Clockwise & counter clockwise movement
Fwd 150° - Rev 30° Fwd 180° - Rev 30°
Fwd 210° - Rev 30° Fwd 250° - Rev 30°
For all continuous rotation files



Nine Editable Memory Programs (M1 to M9)

- RPM Maximum Speed 1000 rpm
- Forward Direction of Rotation
- 360° 360° Clockwise Rotational Movement
- Reverse Direction of Rotation
- Editable Speed / Rpm Setting
- Torque Auto Reverse
- Editable Torque setting



One Piece Contra Angle Head

- 360° Adjustable Handpiece
- Detachable Head
- Miniature Head
- Autoclavable



E-connect S

Endomotor With Built-in Apex Locator



Torque Calibration
For precise torque output




OLED Screen
Digital OLED Display



Ambidextrous Display
Screen can be rotated as per Left or right hand use



Powerful Battery
Heavy duty 1500mAh Lithium battery



M0 - Stand Alone Apex Locator Mode



Apex Locator Bar



Zoom-in & Flash Bar
(when file goes beyond 0.5mm)

APEX LOCATOR

- **Apical control with rotary instruments**
Position of file is displayed on screen and the motor rotates in reverse direction when file reaches pre-set apical working length, to prevent over instrumentation
- **Manual canal measurement**
Manual file clip can be used to measure canals precisely with hand files
- **Multi frequency technology with high accuracy & precision**
- **Apex locator flash bar with zoom-in function**

Working Modes:

FWD - Clockwise Continuous Rotation

REV - Anti-clockwise Continuous Rotation

REC - Reciprocation Movement

- First-in-class cordless open motor with adjustable angles of rotation ranging from 30° to 370° for reciprocation as per user's need

- First-in-class feature of speed adjustment in reciprocation mode

ATC - Adaptive Torque Control


ATC manages the file stress in a smart way. When torque limit is reached, instead of going into reverse mode, motor reciprocates till the stress is reduced. Once back to below pre-set torque level, file will again rotate in continuous rotation

Integrated mode features:

- **Auto Start & Stop** - Motor starts working when the file enters inside root canal, and stops when withdrawn
- **Apical Reverse** - When the file reaches pre-set apical reverse point, user can define whether the file starts moving in reverse direction
- **Apical Slow Down** - To preserve the cemento-dentinal junction and avoid file separation, speed of motor slows down when file is in apical region

Specifications:

Rotational Speed Range : 120-1000 rpm
Torque Range : 0.5-4.0 N.cm
Memory Programs : 11 (M0 to M10)
Battery Capacity : 1500 mAh



E-connect S+

Endomotor With Built-in Apex Locator



360° Rotatable Contra Angle

Upgraded contra angle
One-piece body design



TFT Colorful Display

With excellent color saturation, tint
reducing power and higher contrast



Torque Calibration

For precise torque output



Seamless Shell

Ergonomic design



Advanced Brushless Motor Less Reciprocal Vibration

With reliable quality, bring you precise
operation and higher efficiency



Visualized Apical Display

A simulated image of the apex will appear on the
screen when using the apex locator function, which
vividly shows the apex position and is convenient for
dentists' judgment



Ambidextrous Display

Screen can be rotated, convenient for lefties



Adjustable Backlight Brightness

The brightness of backlight is adjustable
which can save power and meet the
photosensitive needs of different dentists




Built-in Apex Locator - Sensitivity Adjustable


The sensitivity of apex locator is adjustable

Customizable Built-in Endo File System


User could input endo file information



Built-in File System



Long Press **S** to enter
the Mode Setting interface,
press "New File" to add a new file system





Long press **S** to modify the File name






Taper Indicator	Torque Limit	Speed Set	Operation Mode	Size	Taper
	2.5 N·cm	250 rpm	FWD	19	2

Press to set Taper, Size, Operation Mode, Speed Set, Torque Limit, Taper Indicator



"New file" mode can be chosen to use



Can also continue to add new file parameters in this mode

Working Modes:

- **Fwd** - Forward Rotation
- **REC** - Reciprocation Mode
 - Reciprocating angle range 30°-370°
 - The reciprocating angle can be adjusted every 5 degrees
- **Rev** - Reverse Rotation
- **ATC** - Adaptive Torque Control Mode
 - After reaching the preset torque, motor will reciprocate according to the set angle to improve the efficiency of root canal enlargement



Specifications:

- Speed Range**
50 rpm -1500 rpm
- Torque Range**
0.5-5.0 N·cm
- Reciprocating Angle**
Every 5 degrees
- Battery Capacity**
1900 mAh



E-value

Cordless Endo Motor

Specifications:

- Speed Range : 120-1000 rpm
- Torque Range : R.L-4.0 N.cm
- Memory Programs : 10 (M1 to M10)
- Battery Capacity : 1900 mAh



OLED Color Screen
More comfortable visual experience



Built-in File System
Switch file systems quickly



360° Rotatable
360° Rotatable Contra Angle



Brushless Motor
Reliable, Durable, Efficient



Long Battery Life
Heavy capacity
1900mAh lithium battery



Mini Contra Angle
Advanced 1:1 contra-angle
One piece design




New contra-angle head



Normal head & neck


Built-in File System



Note: the doctor could adjust the speed, torque limit and reciprocating angle according to the needs of the case.

Working Modes:

- **Fwd** - Forward Rotation
- **Rev** - Reverse Rotation
- **REC** - Reciprocation Mode
Doctor could set the reciprocating angle and speed level according to actual needs
- Speed Level range from Level 1 to Level 5
- Angles range from 30° to 370°
- **ATC** - Adaptive Torque Control Mode
After reaching the preset torque limit, motor will reciprocate according to the set angle to improve the efficiency of root canal enlargement



Ultra X


Cordless Ultrasonic Activator

Available Tips

Size	Length(mm)
#20 2%	18 & 21
#25 2%	18 & 21



- Autoclavable Tips
- Smooth Surface Finish
- Flexible Tips
*Only X Blue & Silver are Flexible Tips
- Multiple Tip Options




In Action

Cavitation Effect

It is the process of formation of millions of gas bubbles in the irrigation solution. These gas bubbles grow in size, oscillate in different directions and explode into the liquid, producing powerful cleaning action

Acoustic Streaming

It is a unidirectional flow of liquid produced by ultrasonic waves. Movement of tip produces special effect called acoustic turbulence which accelerates the movement in liquid inside the root canal




Indication

Smear Layer Removal

Effective in removal of smear layer and unplugging of dentinal tubules

Dentin Chips Removal

Ultrasonic agitation forms suspension of residual dentinal chips and tissue remnants with irrigant and facilitates their easy removal




Indication

Biofilm Disruption

Ultrasonic cavitation disrupts the biofilm and helps to reduce bacteria load in the canal

Cleaning of Complex Anatomy

Ultrasonic activation helps irrigants to flush into undetected or un-instrumented canal anatomy like isthmus, fins, lateral canals, apical ramification etc. Also helps to remove vapour lock



Protocol

- Prepare the root canal up to required size with the help of shaping and finishing instruments
- Choose the desired disinfectant and flushing agents for cleaning, debridement and disinfection of canals
- Fill the chosen irrigation solution into the canal, never run Ultra X without irrigating solution
- Select the Ultra X tip of desired size and length which passively fits into canal, always keep 2 - 3mm short of apex
- Activate solution using short vertical up & down strokes for 30 seconds. Be careful to not breach apical 2mm
- Use the device gently and never apply pressure or wedge the tip inside the canal. If any points of resistance, then retrieve the ultrasonic tip
- Replenish the irrigating solution, use suction to remove the loose debris. Repeat the cycle 4 to 5 times
- Completely dry the canal and fill it with choice of obturation material
- It is always recommended to isolate the tooth using rubber dam and use personal protective barriers like head gear, mask, eyeglasses etc

Key Features

- Ultrasonic Activation Device
- 45 kHz High Frequency
- Application Based Multiple Tips
- 1500 mAh Powerful Battery

Key Features

- Two Power Modes with LED Indicator
- Contra-angle Design
- Light Weight & Ergonomic
- Cordless Irrigation Device

Battery


- 1500 mAh Lithium Battery
- Battery Level Indicator
- 4.5 Hrs of Continuous Operation
- Option of Direct Plug & Play



UltraMint & UltraMint Pro

All-round Ultrasonic Scaler


Scaling | Periodontic | Endodontic



Soft Perio Power for Ultimate Comfort

Autoclavable plastic box - tips*6, universal wrench*1, EMS type handpiece*1

Autoclavable plastic box - tips*6, universal wrench*1, Satelec type handpiece*1



*Applicable for all tips from Eighteenth


Smart touch system-ultrasensitive & waterproof



Intelligent self-cleaning-makes disinfection easy and avoids cross-infection





Two ways of water supply- Changing supply methods is easier than ever before



900ml stand-alone bottle with clear 20ml scale mark makes it convenient for formulating liquid.




Detachable handpiece with LED- autoclavable, ergonomic, brighter visibility



Detachable and autoclavable handpiece holder, resists cross-infection



Water control directly from handpiece



Ergonomic foot pedal



• EMS type

• Satelec type

One scaler, 2 solutions

Wide range of frequency: 25kHz-42kHz

- We are proud to give you two solutions to perform at your best- both handpieces (EMS type and Satelec type) provided by Eighteeth are compatible with UltraMint & UltraMint Pro.* Scaling can be as convenient as you can imagine.

*Satelec and EMS are trademarks of Acteor® and EMS®, respectively.

AI Engine – Scaling has never been so intelligent!



Customized treatment based on tooth surface conditions



Smart touch control system- ultrasensitive & waterproof



Elliptical stroke patterns with polishing effect




Intelligent degenerative feedback



Digital controlled automatic frequency tracking
Ensures machine works on optimum frequency



Sine wave drive
Avoids simulation caused by instantaneous current output; ensures comfortable treatment



Tips: G1 G2 G4 P1 P3 E4
universal wrench, EMS type handpiece



Tips: E3-SD E5-S E7-S P4-S G1-S G2-S
universal wrench, Satelec type handpiece

Periodontic



P1	P1-S	P3	P3-S
Removal of subgingival calculus and periodontal irrigation		Removal of calculus deep in periodontal pocket and its irrigation	
P4	P4-S	P52	P52-S
Removal of calculus in shallow pockets		Supragingival scaling	
P50L	P50L-S	P50R	P50R-S
Used for subgingival calculus probing, subgingival proximal surface scaling and root planing			
P53L	P53L-S	P53R	P53R-S
Used for subgingival calculus probing, subgingival proximal surface scaling and root planing			
P59	P59-S	P56	P56-S
Used for subgingival calculus probing and scaling		Removal of significant supragingival calculus and plaque	
P3D	P3-SD	P4D	P4-SD
Used to level off the surface of endo after the periodontal flap surgery		For locating root canal and removal of calcification from one-third of crown	

Endodontic



E1	E1-S	E3	E3-S
120° angle holder for 0.7mm-0.8mm endo file for root canal cleaning		For lateral condensation of gutta percha by heating effect	
E4	E4-S	E5	E5-S
Removal of foreign filling material during retreatment			
E6	E6-S	E7	E7-S
Removal of foreign material in cavity to access canal orifice		Removal of broken instrument	
E8	E8-S	E9	E9-S
Flexible and used to brush the canal walls during retreatment		Removal of broken instrument	
E10	E10-S	E3D	E3-SD
Used for root planing in the retrograde preparation of root canal		Enlargement of the canal walls to remove the broken instrument	
E30	E30-S	E30-SD	
Removal of deep foreign particles in canal			

Endodontic



E31	E31-S	E31-SD
Removal of foreign particles from one-third of the canal		
E4D	E4-SD	E5D
Removal of hard material in canal during retreatment		Removal of hard material in canal during retreatment
E6D	E6-SD	E7D
Removal of foreign material in cavity to access canal orifice		Removal of calcification and filling material in cavity to access canal orifice during retreatment
E8D	E8-SD	E9D
Used for root canal retrogression, efficient root apical polishing		Removal of broken instrument
TC1	TC1-S	TC1D
Removal of cervical caries and caries on the occlusal surface		

Scaling



G1	G1-S	G2	G2-S
Removal of supragingival calculus		Removal of significant supragingival calculus	
G4	G4-S	G5	G5-S
Used for supragingival and interdental space scaling		Removal of marginal supragingival calculus and plaque	


*Satelec and EMS are trademarks of Acteon® and EMS®, respectively.

*Satelec and EMS are trademarks of Acteon® and EMS®, respectively.

Fast-Pack pro & Fast-Fill

3D Obturation Devices






Fast-Pack pro


Obturation Pack Handpiece




Temperature Range

-  - Wide temperature range of 90° to 300°
- Five memory programs to customize temperature and cut-off time as per obturation material used or any other intended use of the device


Internally Heated Tip

-  - Only the working end (4 - 6 mm) of the tip is heated - avoids need for insulating sleeves
- The non-working end is insulated and stays cool avoiding soft tissue burn


Class-leading Controlled Heating

-  - Fast heating - Quick heating prevents PDL damage
- Swift cooling - Swift cooling ensure GP does not follow along the heating tip

Experience Ergonomics

-  - 360° start / stop ring button allows holding the device in any angle comfortable for you
- Left - right handed operation
- Comfortable grip, light in weight

Powerful Battery

-  - Highly efficient 2600mAh Lithium battery
- Direct Plug & Play - In case of low battery, use the device directly with adaptor

Specifications

Heating Time	0.2 seconds for 200 °C Less than 1 second for 300 °C
Cooling Time	2 seconds for cooling
Charging Time	2 hours to complete charge
Operating Time	4 hours of continuous operation at 250 °C





Fast-Fill

Obturation Fill Handpiece



Fast-Fill Tips



- Autoclavable silver tips for multiple usages
- Silver tips maintain the thermoplastic temperature of GP and offers smooth flow inside canal
- 360 ° rotating tip for easy filling at all angles
- Extra long needle offers convenience and visibility
- Easy to maintain and replace

Motorized GP Flow



- Provides consistent flow of GP at desirable rate
- Precludes voids and makes obturation homogeneous
- Eliminates hand fatigue that resulted from earlier gun type devices

Comprehensive Display



- Three customizable memory programs - T1 to T3
- Temperature range 100 °C - 200 °C
- GP flow rate - Low, Medium, High
- Reserve GP level indicator

Intuitive Design



- Sleek and ergonomic hand-piece
- Simple push button to dispense GP
- LED indicator reflects when the device is ready for use

Powerful Battery



- Highly efficient 2600mAh Lithium battery
- Direct Plug & Play - In case of low battery, use the device directly with adaptor

Specifications


Heating Time	20 seconds to 200 °C
Charging Time	4 hours to complete charge
Hexagonal Wrench	Use the wrench to bend the tip, rotate the tip and clean the device



23 Gauge Large



25 Gauge Small



Protocol for Use

Master Cone Verification & Resizing

Select a master cone which corresponds to the master apical preparation. Check for tug-back and confirm the same with radiograph. Resize the GP by cutting it to 0.5 to 1mm short of the working length. This will prevent peri-apical extrusion of the GP

Selection of Fast-Pack & Plugger Tip

Select the Fast-Pack tip which reaches 4-5 mm short of the working length and give a snug fit at the tip. Mark this length using the silicon stopper. Similarly, select the hand plugger which reaches 4-5mm short of the working length

Gutta Percha & Sealer

Coat the master cone apically with the sealer and insert it into the canal. Brush the canal wall using strokes of GP and sealer



Down Pack with Fast-Pack

Remove the GP from the orifice/coronal part with the help of heated Fast-Pack tip. Plug the GP cone with pre-fit hand plugger at the orifice. Now take the heated Fast-Pack tip 4-5mm short of the working length. The bite of GP will come out of the canal along with the tip. Heated Fast-Pack tip will transfer heat to the GP and convert it into thermoplastic state. Condense the GP with the help of hand plugger to ensure the complete sealing of accessory canals and apex

Fill with Fast-Fill

Filling the canal using the Fast-Fill device in increments of 4mm till the canal is filled completely while condensing the warm GP with the appropriate plugger after each increment


Gutta Percha Points




Ideal for obturating root canals which have been prepared with GT files, Profile, Quantecs, etc. Greater Taper 0.04 and 0.06 are available in assorted or individual sizes. At the top of points (d1), the diameter corresponds to the ISO size, but then the conicity is greater than ISO points.

	NO.15 NO.20 NO.25 NO.30	○ white ● yellow ● red ● blue	GP 15 20 25 30
	NO.20 NO.25 NO.30 NO.35	● yellow ● red ● blue ● green	GP 20 25 30 35
	NO.20 NO.25 NO.30 NO.35	● yellow ● red ● blue ● green	GP 20 25 30 35
	NO.25 Gutta Percha Points (for Pro Taper Gold) F2	● red	GP 25
	NO.30 Gutta Percha Points (for Pro Taper Gold) F3	● blue	GP 30

Absorbent Paper Points



	NO.15 NO.20 NO.25 NO.30	○ white ● yellow ● red ● blue	PP 15 20 25 30
	NO.20 NO.25 NO.30 NO.35	● yellow ● red ● blue ● green	PP 20 25 30 35
	NO.20 NO.25 NO.30 NO.35	● yellow ● red ● blue ● green	PP 20 25 30 35
	NO.25 Absorbent Paper Points (for Pro Taper Gold) F2	● red	PP 25
	NO.30 Absorbent Paper Points (for Pro Taper Gold) F3	● blue	PP 30



CuringPen

4 LEDs Curing Light



Technical parameter:

Dimension: 24*21*230(mm)

Net Weight: 109g

Lens size: 11mm

Power supply:

>1200 times, 3 Sec

Maximum light intensity




Smart controlling chip
keeps light output consistent

Light intensity:

1000mw/cm², 1500mw/cm², 2300mw/cm²


600 mw/cm²(detect mode)



Wavelength:385-515nm
3 Blue, 1 Violet, 10w LED



360° rotatable head




Intensity digital display




Multiple Curing Options


Ramp mode



Pulse mode



Caries detection





CuringPen-E

3 LEDs Curing Light


Technical Specifications:

Dimension	: 19 x 16.5 x 7.5 cm
Weight	: 120 g
Lens Size	: 8 mm
Power Supply	: >1200 times, 3000 mw/cm ² , 3 Sec
Maximum Light Intensity	: 3000 mw/cm ²
Spectral Range	: 380-515 nm



360° Rotatable Head

Easy to illuminate various treatment areas



Small Light Head

Easy to access to posterior molars which offers better patient experience




Special Lens Design

Unique lens design forms a concentrated beam for uniform curing



Ultra Power

- Light intensity now up to 3000mw/cm²
- High light intensity results in deeper cure depths which could curing resins more efficiently and with short time



Broad Spectral Range 380 nm-515 nm

- Why broad spectral range?
- Most of resins use camphorquinone (CQ) materials, but some manufacturers choose alternative initiators such as TPO to reduce the yellow impact caused by CQ. TPO can only be activated by the light wave below 410nm. CuringPen-E has a broad spectrum which could effectively cure different resins.

Working Modes

• P1



Normal Mode

Light Intensity : 1200mw/cm²

Time Setting : 5s, 10s, 15s, 20s, 25s, 30s, 35s, 40s

• P2



High Mode

Light Intensity : 2000mw/cm²

Time Setting : 1s, 3s, 5s, 10s

• P3



Ultra Mode

Light Intensity : 3000mw/cm²

Time Setting : 1s, 3s, 5s

• P4




Pulse Mode

Light Intensity : 1200mw/cm²

Time Setting : 5s, 10s, 15s, 20s

Pulse mode



• P5





Ramp Mode

Light Intensity : 1200mw/cm²

Time Setting : 5s, 10s, 15s, 20s

Ramp mode






HyperLight

Portable X-ray Unit




Engineered for accurate images

- Clinical Imagery



- Sharp, high resolution images with a 0.4 focal spot

Source focal spot Source focal spot



Effect of focal spot size



Designed with safety in mind



Built-in leaded double shield to protect the operator from leakage radiation



Locking system to prevent unintentional exposure



Specifications

Model	HyperLight	Weight	1.8 kg
Focal Spot	0.4	Battery	2500 mAh x 4
Tube Voltage	65 kV	Exposure Time Range	0.02-2 Sec
Tube Current	2.5 mA	Line voltage	100V-240V, 50Hz/60Hz




Simplify your digital workflow


- Compact size and lightweight can make you operate it by one hand only
- Intuitive user interface
- Smart angulation system



Maximize your diagnostic potential

Perfect combination





HyperLight-G

Portable X-ray Unit



Receptor

Use with film, digital sensor, phosphor plate



Locking System




Preset Parameter

Divided into adults and children, and can be set and controlled for different teeth



Large Battery Capacity

More than 300 exposures in one full charge



Smart angulation system, faster your workflow



Safe

- Internal x-ray shielding
- Locking system to avoid unintentional exposure
- Backscatter shield for operator and patient protection(*optional)

Portable

- Cordless, handheld and extremely mobile, move freely between rooms, locations and cities

Economical

- Use with film, digital sensor or phosphor plate
- Prolonged working life, at least 80,000 exposures

Specifications

Model	HyperLight-G
Tube Voltage	65 kV
Tube Current	2 mA
Exposure Time Range	0.02-2 Sec
Weight	1.9 kg
Battery	2000mAh X4
Line Voltage	100V-240V, 50Hz/60Hz